











# Agilent Vial Identification Chart

## with Recommended Fill Volumes




### Wide Opening Screw Top Vials (9 mm)      Narrow Opening Screw Top Vials (8 mm)      Polypropylene Vials

									
32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm
1.5 mL	1.2 mL	250 µL	1.3 mL	1.5 mL	250 µL	250 µL	250 µL	250 µL	0.7 mL
5182-0714	5183-2030	5188-6591	5184-3550	5183-4428	5188-2788	9301-0977	5190-2242	5188-5390	5182-0567
5182-0717	5182-0717	5182-0552	5182-0552	5183-4442	5182-0552	5182-0552 Glass Insert	5182-0717	5182-0717 Glass Insert	5182-0552




### Wide Opening Crimp Top Vials (11 mm)

						
32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm
1.5 mL	1.2 mL	0.7 mL	250 µL	250 µL	1.3 mL	1.2 mL
5181-3375	5182-3454	5182-0567	9301-0978	9301-0977	5184-3551	5185-5821
5182-0552	5182-0552	5182-0552	5182-0552	5182-0552	5182-0552	5182-0552






### Inserts for Wide Opening Vials (11 mm & 9 mm)

		
30 mm x 5.6 mm	31 mm x 5.6 mm	31 mm x 5.6 mm
250 µL	250 µL	300 µL
5181-1270	5183-2085	5181-3377





### Inserts for Narrow Opening Vials (8 mm)

		
28 mm x 4.8 mm	31 mm x 4.8 mm	31 mm x 4.8 mm
150 µL	150 µL	150 µL
5183-2088	5183-2089	5183-2090



### Wide Opening Snap Top Vials

				
32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm	32 mm x 12 mm
1.7 mL	1.4 mL	0.7 mL	0.3 mL	0.3 mL
5182-0544	5183-4510	5182-0567	9301-0978	9301-0977
5182-3458	5182-3458	5182-3458	5182-3458	5182-3458




### Headspace Vials

			
75 mm x 22 mm	45 mm x 22 mm	75 mm x 22 mm	45 mm x 22 mm
20 mL	10 mL	20 mL	10 mL
5182-0837	5182-0838	5183-4474	5183-4475
5183-4477	5183-4477	5183-4477	5183-4477

### Mini Vials

	
500 µL	700 µL
5180-0806	5180-0805
5180-0842	5180-0842

### Specialty Vials

		
37 mm x 22 mm	45 mm x 15 mm	37 mm x 22 mm
6 mL	4 mL	6 mL
9301-1419	5183-4448	9301-1377
9301-1425	5183-4464	9301-1379